



### Deluxe Automatic Overbed Table - Model OBT500

The Deluxe Automatic Overbed Table features durable, heavy duty construction and an attractive functional design.

#### Features:

- Adjustable from 30" to 45", in one-inch increments, with the use of a counterbalanced spring height adjustment
- Heavy-duty U-base with 28" clearance and 2" twin wheel casters
- Extra wide high pressure laminate top with T-mold edge measures 18" x 32"
- Capacity - 100 lbs. evenly distributed
- Base available in chrome, or optional powder coat finishes
- Specify Joerns laminate finish for top

#### Accessories:

Double Mirror Vanity

#### Joerns OBT500

#### Laminate Top Finishes:

- Beigewood
- Natural Maple
- Provincial Cherry
- Quartered Oak
- Select Cherry
- Valley Oak
- Walnut
- Wild Cherry



### Automatic Overbed Table - Model OBT200

The Automatic Overbed Table offers an adjustable tabletop that can be raised and lowered with the touch of a finger.

#### Features:

- Adjustable from 29½" to 45", using a mechanism that automatically locks in place at the desired height
- High pressure laminate top with T-mold edge measures 15" x 30"
- U-base
- Deluxe 1½" heavy duty casters standard
- Capacity - 50 lbs. evenly distributed
- Base available in Brown or Soft Tone
- Specify Joerns laminate finish for top

#### Accessories:

Double Mirror Vanity

#### Joerns OBT200

#### Laminate Top Finishes:

- Beigewood
- Natural Maple
- Provincial Cherry
- Quartered Oak
- Select Cherry
- Valley Oak
- Walnut
- Wild Cherry



### Automatic Low Overbed Table - Model OBT250

The Automatic Low Overbed Table is designed specifically for use with low beds, and features an adjustable table top that can be quickly and easily adjusted.

#### Features:

- Infinitely adjustable from 20¾" to 27¾"
- High pressure laminate top with T-mold edge measures 15" x 30"
- U-base with deluxe 1½" heavy duty casters standard
- Capacity - 50 lbs. evenly distributed
- Base available in Brown or Soft Tone
- Specify Joerns laminate finish for top

#### Accessories:

Double Mirror Vanity

#### Joerns OBT250

#### Laminate Top Finishes:

- Beigewood
- Natural Maple
- Provincial Cherry
- Quartered Oak
- Select Cherry
- Valley Oak
- Walnut
- Wild Cherry



**Deluxe Automatic Overbed Table - Model OBT300**

The Deluxe Automatic Overbed Table features durable, heavy duty construction and an attractive functional design.

**Features:**

- Infinitely adjustable from 28" to 45"
- H-base with 2" oversized casters
- Designer veneer top measures 15" x 30"
- Top features a ¼" spill protection lip
- Capacity - 50 lbs. evenly distributed
- Base available in Soft Tone or optional powdercoat finishes
- Specify laminate finish for top

**Accessories:**

Double Mirror Vanity

**Joerns OBT300**

**Laminate Top Finishes:**

- Fusion Maple
- Grey Nebula
- Gunstock Walnut
- Honduras Mahogany
- Natural Pear
- Quartered Oak
- Wild Cherry



**Automatic Overbed Table - Model OBT100**

This Automatic Overbed Table features a basic design to meet the requirements of healthcare facilities.

**Features:**

- Adjustable from 28" to 42"
- Low-pressure laminate top with T-mold edge measures 15" x 30"
- H-base with 1½" twin wheel casters
- Capacity - 50 lbs. evenly distributed
- Base available in zinc only
- Top available in walnut only