



MANGAR Leglifter

The inflatable Leglifter raises, lifts and transfers the users' legs into bed.

EASY TO FIT

- Fits most beds including profiling beds
- Straight forward to fit or remove, and requires no permanent installation
- Attaches to the bed underneath the mattress
- One airline connects the Leglifter to the Airflo Compressor and inflation is activated by pushing one button

SIMPLE TO MAINTAIN

- For a routine clean, sponge clean with a detergent solution and dry thoroughly
- The covers can be machine washed
- Spare covers can be purchased

- The materials are designed to exceptionally high standards with extraordinary durability, strength & performance in the toughest of environments

SAFE TO USE

- Raises, lifts and transfers legs into the bed
- Leg break on profiling beds can still be used once Leglifter is fitted
- Aids independence

IMPROVES WELL-BEING

- Minimizes moving and handling so reducing the risk of musculoskeletal injury for the caregiver
- Fits discreetly against the side of the bed
- No permanent installation required

TECHNICAL SPECIFICATIONS

Max user weight	297 lbs
-----------------	---------

Min height	Flat
------------	------

Total length	49"
--------------	-----

Width	22.5"
-------	-------

Total product weight	3.25 lbs
----------------------	----------


Airflo Compressor	12
-------------------	----

Lifting performance per charge	Up to 16 lifts
--------------------------------	----------------



Features

1. Fitting Straps.
2. Air Distribution Hose.
3. Made from strong polyurethane coated nylon material.
4. Seams joined by high frequency welding.
5. Hook of Leglifter straps.



Product code: MPCA140500



Joerns.com

General 800.826.0270

Canada 866.546.1151

6240088 RevA • 000576

US Partner with

